

NITRILE

E3NA Accelerator Free Powder Free Nitrile Examination Gloves

Non - Sterile

Non Sterile • Powder Free • Ambidextrous • Single Use • Beaded Cuff • Fingertip textured • Single Chlorination • Synthetic Rubber Nitrile

Facts & Features:

- Compound to meet our cleanliness criteria
- Cuff strengthened by a rolled cuff
- Easy donning and removal
- Excellent tactile sensitivity (Micro textured fingertips)
- Excellent viral barrier properties
- Free of Accelerator
- Free of Latex
- Free of Powder
- Free of Protein
- Free of thiuram, plasticizers, animal tissue and other biological substances
- Good sensitivity and feel
- Hot water leached
- Low risk personal protective use
- Low dermatitis
- Online Chlorination
- Safe for food handling

Dental
Food
Healthcare
Laboratory
Hairdressing
Personal Protection

Colours Availability:



White



Blue

Product Specifications:

Size	Article No.	Palm Width (in mm)	Force @ Break (Newton)	
			Throughout Shelf Life	After Challenge
XS	E3NA01	75±5	≥6N	≥6N
S	E3NA02	85±5		
M	E3NA03	95±5		
L	E3NA04	106±5		
XL	E3NA05	116±5		
Cuff Thickness		≥0.05	Available Length	≥240mm / ≥295mm
Palm Thickness		≥0.06	Shelf Life	3 Years
Finger Thickness		≥0.08	Weight M size	4.0±0.2g / 5.0±0.2g

Conformity standards/ Regulation:

1. Class I Medical Device as per MDD 93/42/EEC (Conform to EN 455 part 1, 2, 3 and 4 and ISO 10993 Biological Evaluation)	CE
2. PPE Directive 89/68/EEC – CAT III, (Conform to EN 374, EN 388 and EN 420)	CE 0321, Biohazard, Sharp, Flammable
3. Product compliance with Regulation (EC) no 1935/2004 & meeting demand of Commission Regulation (EU) no 10/2011	Food Contact
4. Produced under ISO 13485:2003 and EN ISO 13485:2012 Quality Management	BSI, UKAS
5. Produced under ISO 9001:2008 Quality Management	BSI, IAF, ANAB

Email us for further information:

EN Medipart (Asia Pacific) Sdn. Bhd. • E-Mail: info@en-medipart.com.my • www.en-medipart.com.my